



SHELDON STORAGE CENTERS

Sheldon Storage Center

is the latest product innovation from Sheldon Labs. Offered in a variety of sizes and tote tray options, these cabinets are the perfect storage solution for your classroom.

Features

- Durable SEFA laboratory-grade, wood cabinet construction including 5-knuckle hinges.
- Tote trays are shatterproof storage that combine strength, durability, and design.
- Cabinets are available with or without doors and offered as stationary or mobile.



sheldonlabs.com

STORAGE CENTER CABINETS



TALL CABINET

Available in two sizes:
83"H x 28"W or 83"H x 41.5"W



BASE CABINET

Multiple Sizes Available
Heights: 36" or 33"
Widths: 41.5", 28", or 14"



MOBILE CABINET

Two Styles Available
Mobile Cubicle (pictured):
24.5"H x 15"W x 18.5"D

Mobile Base Cabinet:
33"H x 29"W

All cabinets available with and without doors.

STORAGE CENTER TRAYS



SMALL TRAY



MEDIUM TRAY



LARGE TRAY



GLIDE & TILT TRAY

Trays are available in 3 sizes: small, medium, and large. The glide and tilt trays can be opened in an angled, open position making it easy to access the contents. The runner holds the tote in place and prevents it from falling out of the unit. When necessary, the tote is easily removed from the unit to transport elsewhere. Tote includes a window to insert your own label.

Sheldon Laboratory Systems

102 Kirk Street | Crystal Springs, MS 39095 | 800.531.7604