

A. **Sheldon TM426122-200 Privacy Panel Storage Cabinet:**
42" L. X 22" W. X 60-1/2" H.

1. The Privacy Panel Storage Cabinet is constructed of 7-ply hardwood plywood with double doors and finished back. The interior of the unit is furnished with three (3) horizontal rails that support up to twenty one (21) Axis table privacy panels. Unit includes one (1) adjustable shelf and twelve (12) tote trays for apparatus storage. Four (4) locking swivel castors are shipped separate and field installed. The unit is available in Oak or Maple, as specified.
2. Standard accessories included with each unit:
Twelve (12) #1700-0112 Tote trays.

Optional accessory available at additional cost.

#PPA0001K Privacy Panel Kit. (three (3) panels per kit; one (1) kit recommended for each four (4) student Axis Table).