

A. **Sheldon 20620 Instructor Demonstration Center:**
60" L. X 31" W. X 36" H. Sink on Left End.

TP 12.14.17

1. Demonstration Center: Includes one (1) S183528-101 sink cabinet, one (1) B183528-571 five-drawer cabinet, one (1) kneespace, and a finished back. The Center also includes a black Shelresin top with epoxy resin drop-in sink on left end. **Internal service piping and wiring not included.**
2. Standard services and accessories included:
 - a. One (1) 80020-CV Unicast combination cold water and double gas fixture with check valve.
One (1) D05 Shelresin 14" x 10" x 6" I.D. drop-in sink.
One (1) sink outlet, strainer and stopper.
One (1) 85101 GFI duplex electrical outlet.
One (1) 86380 upright rod assembly with two (2) 86320 rod bases.

A. **Sheldon 20625 Instructor Demonstration Center:**
60" L. X 31" W. X 36" H. Sink on Right End.

TP 12.14.17

1. Demonstration Center: Includes one (1) S183528-100 sink cabinet, one (1) B183528-571 five-drawer cabinet, one (1) kneespace, and a finished back. The Center also includes a black Shelresin top with epoxy resin drop-in sink on right end. **Internal service piping and wiring not included.**
2. Standard services and accessories included:
 - a. One (1) 80020-CV Unicast combination cold water and double gas fixture with check valve.
One (1) D05 Shelresin 14" x 10" x 6" I.D. drop-in sink.
One (1) sink outlet, strainer and stopper.
One (1) 85101 GFI duplex electrical outlet.
One (1) 86380 upright rod assembly with two (2) 86320 rod bases.

A. **Sheldon 20630 Instructor Desk Extension:**
50-1/2" L. X 31" W. X 30" H. Right Hand.

TP 12.14.17

1. Desk Extension: Includes one (1) B182928-500 three-drawer cabinet, one (1) knee-space, and a finished back. The cabinet is furnished with two (2) small drawers and one (1) file drawer. The assembly also includes a black Shelresin top. Cabinet is located on right end.

A. **Sheldon 20635 Instructor Desk Extension:**
50-1/2" L. X 31" W. X 30" H. Left Hand.

TP 12.14.17

1. Desk Extension: Includes one (1) B182928-500 three-drawer cabinet, one (1) knee space, and a finished back. The cabinet is furnished with two (2) small drawers and one (1) file drawer. The assembly also includes a black Shelresin top. Cabinet is located on left end.

A. **Sheldon 20650 Instructor Demonstration Center:**
96" L. X 31" W. X 36" H. Sink on Left End.

TP 12.14.17

1. Demonstration Center: Includes one (1) S353528-200 sink cabinet, one (1) B353528-570 seven-drawer cabinet, one (1) 24" kneespace, and a finished back. The Center also includes a black Shelresin top with epoxy resin drop-in sink on left end. **Internal service piping and wiring not included.**
2. Standard services and accessories included:
 - a. One (1) 80020-CV Unicast combination cold water and double gas fixture with check valves.
One (1) D20 Shelresin 16" x 16" x 7-1/2" I.D. drop-in sink.
One (1) sink outlet, strainer and stopper.
One (1) 85101 GFI duplex electrical outlet.
One (1) 86380 upright rod assembly with two (2) 86320 rod bases.

A. **Sheldon 20655 Instructor Demonstration Center:**
96" L. X 31" W. X 36" H. Sink on Right End.

TP 12.14.17

1. Demonstration Center: Includes one (1) S353528-200 sink cabinet, one (1) B353528-570 seven-drawer cabinet, one (1) 24" kneespace, and a finished back. The Center also includes a black Shelresin top with epoxy resin drop-in sink on right end. **Internal service piping and wiring not included.**
2. Standard services and accessories included:
 - a. One (1) 80020-CV Unicast combination cold water and double gas fixture with check valves.
One (1) D20 Shelresin 16" x 16" x 7-1/2" I.D. drop-in sink.
One (1) sink outlet, strainer and stopper.
One (1) 85101 GFI duplex electrical outlet.
One (1) 86380 upright rod assembly with two (2) 86320 rod bases.

A. **Sheldon 20672 ADA Instructor Demonstration Center:** TP 12.14.17
72" L. X 31" W. X 33" H. Sink on Left End.

1. Demonstration Center: Includes one (1) S183228-101 sink cabinet, one (1) B183228-540 four-drawer cabinet, one ADA (1) kneespace, and a finished back. The Center also includes a black Shelresin top with epoxy resin ADA drop-in sink on left end. **Internal service piping and wiring not included.**
2. Standard services and accessories included:
 - a. One (1) 80020-WB-CV Unicast combination cw/gas fixture with wrist blade and check valves.
One (1) A05 Shelresin 14" x 10" x 5" I.D. drop-in sink.
One (1) sink outlet, strainer and stopper.
One (1) 85101 GFI duplex electrical outlet.
One (1) 86380 upright rod assembly with two (2) 86320 rod bases.

A. **Sheldon 20674 ADA Instructor Demonstration Center:** TP 12.14.17
72" L. X 31" W. X 33" H. Sink on Right End.

1. Demonstration Center: Includes one (1) S183228-100 sink cabinet one (1) B183228-540 four-drawer cabinet, one (1) ADA kneespace, and a finished back. The Center also includes a black Shelresin top with epoxy resin ADA drop-in sink on right end. **Internal service piping and wiring not included.**
2. Standard services and accessories included:
 - a. One (1) 80020-WB-CV Unicast combination cw/gas fixture with wrist blade and check valves.
One (1) A05 Shelresin 14" x 10" x 5" I.D. drop-in sink.
One (1) sink outlet, strainer and stopper.
One (1) 85101 GFI duplex electrical outlet.
One (1) 86380 upright rod assembly with two (2) 86320 rod bases.

A. **Sheldon 20696 ADA Instructor Demonstration Center:** TP 12.14.17
96" l. x 31" w. X 33" H. Sink on Left End.

1. Demonstration Center: Includes one (1) S343228-000 ADA sink cabinet with removable panel, one (1) B303228-560 six-drawer cabinet, one (1) ADA kneespace, and a finished back. The Center also includes a black Shelresin top with epoxy resin ADA drop-in sink on left end. **Internal service piping and wiring not included.**
2. Standard services and accessories included:
 - a. One (1) 80020-WB-CV Unicast combination cw/gas fixture with wrist blade handle and check valves.
One (1) A25 Shelresin 18" x 15" x 5" I.D. drop-in sink.
One (1) sink outlet, strainer and stopper.
One (1) 85101 GFI duplex electrical outlet.
One (1) 86380 upright rod assembly with two (2) 86320 rod bases.

A. **Sheldon 20698 ADA Instructor Demonstration Center:**
96" l. x 31" w. X 33" H. Sink on Right End.

TP 12.14.17

1. Demonstration Center: Includes one (1) S343228-000 ADA sink cabinet with removable panel, one (1) B303228-560 six-drawer cabinet, one (1) ADA kneespace, and a finished back. The Center also includes a black Shelresin top with epoxy resin ADA drop-in sink on right end. **Internal service piping and wiring not included.**
2. Standard services and accessories included:
 - a. One (1) 80020-WB-CV Unicast combination cw/gas fixture with wrist blade handle and check valves.
One (1) A25 Shelresin 18" x 15" x 5" I.D. drop-in sink.
One (1) sink outlet, strainer and stopper.
One (1) 85101 GFI duplex electrical outlet.
One (1) 86380 upright rod assembly with two (2) 86320 rod bases.