

TE II





The **TE II Student Workstation** combines Sheldon's work and standing areas, allowing students to safely collaborate with each other. Each student has access to more counter space and elbow room with TE II. The workstation is available in many different configurations and sizes to fit your needs.

sales@sheldonlabs.com



800.531.7604

TE II

TE II TABLE MATRIX SHELDON LABORATORY SYSTEMS	MODEL #	REVIT FAMILY	TOP SIZE	HEIGHT	SINK	FIXTURES	FEATURES	STUDENT CAPACITY
	14001	TE II Split Top Full Trough	56" x 56" (2 equal tops with trough between tops)	Adjustable at install 33"–37" H	Trough Sink	Unimix or Unicast	2 outlets-perpendicular to trough	4
	14020	TE II Split Top Full Trough ADA	57" x 61" x 74" (lower ADA top is larger than upper top)	Install ADA height 34" and 32" H	Trough Sink	Unimix or Unicast	2 outlets-perpendicular to trough	4
	14123	TE II Flat Top	56" x 56"	Adjustable at install 33"–37" H	No	No	2 outlets	4
	14041	TE II Flat Top	56" x 67"	Adjustable at install 33"–37" H	#DRS12 Round Resin Sink	Unimix or Unicast	2 outlets	4
	14233	TE II Flat Top SM	48" x 48"	Adjustable at install 33"–37" H	No	No	2 outlets	4
	14241	TE III Flat Top SM	48" x 48"	Adjustable at install 33"–37" H	#DRS12 Round Resin Sink	Unimix or Unicast	2 outlets	4
	14211	TE III Flat Top SM Trough	48" x 48"	Adjustable at install 33"–37" H	Trough Sink #DT36	Unimix or Unicast	2 outlets	4
	14127 ADA	TE III Split Top ADA	74" x 67"	Install ADA height 34" and 32" H	No	No	2 outlets	4
	14001.5	TE III Half with Full Trough	67" x 52"	Height to be determined by perimeter casework	Trough Sink	Unimix or Unicast	2 outlets	2

Rod Bases - All TE II models (except Model # 4001.5) include 4 rod bases in table. Model #14001.5 includes 2 rod bases.

TE II Annex Assembly 141212 – Assembly used to connect the TE II student table to counter top work surface.

sales@sheldonlabs.com



800.531.7604