

A. **Sheldon Slimline Axis Student Table:**

Fixed height, 7ft. Work Surface.

1. The moveable Student Laboratory table shall be multi-purpose and must be able to seat a minimum of four (4) students. The table has a fixed height of 36".
2. The 1" thick Shelresin work surface is curvilinear shaped and measures approximately 84" x 44". The unique contoured shape of the top is provided to eliminate sharp corners from projecting into the aisles.
3. The work surface is supported by a 3/4" sub top, and a pedestal base consisting of hardwood vertical ends and a horizontal support. Vertical ends are attached to two (2) 28" wide steel pedestal feet, mounted on adjustable glides.
4. Color of work surface, sub-top, pedestal base, and feet to be black.